









## THOROUGH EXAMINATION

- High quality optical system offers wide field of view and qualitative images.
- 5 drum steps magnifications allow clear observation of finest details.
- Light tilting feature combined with numerous slit adjustments and filter options provide flexibility in viewing different structures of the eye.
- Integrated yellow filter and background illumination provide optimal examination conditions.

## CUTTING EDGE LED ILLUMINATION

- Brighter and more homogeneous illumination.
- Very comfortable examination for the patient with no heat emission.
- Long lifetime with substantial maintenance cost reduction.

## IMAGE AND VIDEO CAPTURING (OPTION)

• Combined with the DS 550 and Anaeyes software, possibility of capturing images and video sequences.

MAGNIFICATION		SLIT	
Туре	Galilean magnification changer	Slit width	Continuous from 0 mm to 12 mm
Standard Eyepiece	12.5x	Slit length	Continuous from 0 mm to 12mm
Selection	Rotation drum 5 positions	Slit angle	0 to 180° with horizontal scanning capability
Field of view	6x (37 mm)	Tilt	5, 10, 15, 20°
	10x (24 mm)	ILLUMINATION LAMP	LED
	16x (14 mm)		LED
	25x ( 8 mm)	FILTERS	Blue, Red-Free, Neutral (normal use), Red and Yellow
	40x (5.2 mm)		
BASE		DIMENSIONS & WEIGHT	
Longitudinal movement	113mm		255(w) X 320 (D) x 665 (H) mm - 8kg
Lateral movement	107mm	Power-supply	100-120V, 220-240V AC
Vertical movement		Consumption	25 VA
Chinrest (vertical)	66mm	$\mathbf{C}\mathbf{\epsilon}$ Conformity marking	
Options	Polarised filter Hruby Lens Applanation tonometer	DS 550 Digital Camera	1/1.8 CCD type Effective pixel number : 2 MP

As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.



Essilor Instruments USA 8600 W. Catalpa Avenue, Suite 703 Chicago, IL 60656 Phone: 855.393.4647 Email: info@essilorinstrumentsusa.com www.EssilorInstrumentsUSA.com Your local distributor :