



DIGITAL SLIT LAMP

















THOROUGH EXAMINATION

- High-quality optical system offers a wide field of view and qualitative images
- 5 drum steps magnifications allow clear observations of finest details
- Numerous slit adjustments and filter options provide flexibility in viewing different structures of the eye
- Integrated yellow filter and background illumination provide optimal examination condition

CUTTING-EDGE LED ILLUMINATION

- Brighter and more homogeneous illumination
- Very comfortable examination for the patient with no heat emission
- Long lifetime with substantial maintenance cost reduction

IMAGE AND VIDEO CAPTURING (OPTION)

• Combined with the DS 550 and Anaeyes software, possibility of capturing images and video sequences

Magnification		SLIT	
Туре	Galilean magnification changer	Slit width	Continuous from 0mm to 12mm
Standard Eyepiece	12.5x	Slit length	Continuous from 0mm to 12mm
Selection	Rotation drum 5 positions	Slit angle	0 to 180° with horizontal scanning capability
	6x (37 mm)		
Field of view	10x (24 mm)	ILLUMINATION LAMP	155
	16x (14 mm)		LED
	25x (8 mm)	FILTERS	51 5 5 5 1 1 1 1 1 1
	40x (5.2 mm)		Blue, Red-Free, Red and Yellow
Base		DIMENSIONS	
Longitudinal movement	113mm	Table	10.0 (W) x 12.6 (D) x 16.9 (H) in
Lateral movement	 107mm	Power supply	100-120V, 220-240V AC
Vertical movement	30mm	Consumption	25 VA
Chinrest (vertical)	76mm	C € Conformity marking	
OPTIONS	Polarized filters Hruby Lens Applanation tonometer	DS550 Image capture device	1/1.8 CCD type Effective pixel number: 2 MP



Essior instruments USA 08/2018. All rights reserved. As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.